



ClearKlens<sup>TM</sup> IPA

VH1

## Sterile isopropyl alcohol for disinfection in cleanrooms

### Description

ClearKlens IPA is a ready to use disinfectant for the Pharmaceutical cleanroom industry. ClearKlens IPA is based on 70% isopropyl alcohol. ClearKlens IPA is available as 900mL bag in bottle trigger spray, and as 1L bottle or 5L canister. The 900mL bottle has been filtered to 0.2µm and irradiated. The 1L/5L has been filtered to 0.2µm and aseptically filled.

### Key properties

- ClearKlens IPA is a solution of 70% isopropyl alcohol and 30% water for injection (USP quality).
- ClearKlens IPA is a ready to use product for the disinfection of work surfaces and process equipment in cleanrooms
- ClearKlens IPA 900mL trigger spray is an antisuck back- bag in bottle delivery system to eliminate product contamination during use.
- BPR supported pack and formulation
- Double bagged for easy transfer into the cleanroom.

### Benefits

- Ready to use
- 0.2µm filtered
- Gamma irradiated/aseptically filled and double bagged
- Detex labelled for proof of irradiation
- Double bagged for easy transfer into the cleanroom
- Odourless and Non tainting
- Broad spectrum of activity

### Use instructions

ClearKlens IPA 900mL trigger spray is designed for spraying and wiping applications. Hold spray approximately 30-40cm from area to be treated. Spray directly onto surface to ensure complete coverage. Allow sufficient wet contact time before proceeding. Always close nozzle after use.

ClearKlens IPA 1L bottle/5L canister is used for mopping and wiping applications. The 1L bottle can be used with the dedicated touchless dispenser.

For mopping and wiping application dispense ClearKlens IPA directly into bucket. Wet cloth with alcohol and wipe onto surface to ensure complete coverage. Allow sufficient wet contact time before proceeding.





ClearKlens

IPA

VH1

#### Technical data

Appearance: Clear liquid free from insolubles

pH (1% solution at 20°C): 6-9

Density (20°C): 0.870-0.880

Alcohol proportion: 70% +/- 2.5% v/v

*The above data is typical of normal production and should not be taken as a specification.*

#### Safe handling and storage information

Store in original closed containers, away from sunlight and extremes of temperature. Full guidance on the handling and disposal of this product is provided in a separate Safety Data Sheet.

#### Product compatibility

Application solutions, when used as directed, will not affect materials normally encountered in the Pharmaceutical and Cosmetics industries.

#### Test method

Each lot is provided with:

Certificate of Analysis, available also on [www.clearklens.com](http://www.clearklens.com)

Certificate of Sterility

Certificate of Irradiation

#### Microbiological data

ClearKlens IPA has been extensively tested using EN 1276, EN 1650 and EN 13697. These have shown ClearKlens IPA to have a broad spectrum of activity against both Gram-positive/Gram-negative bacteria and yeasts.

Tests according to EN 14476 have shown a biocidal efficacy against enveloped viruses (Modified vaccinia virus Ankara/MVA) like SARS-CoV-2.

#### Available pack sizes

6 x 0.9 L trigger spray

10 x 1 L (bottle, can be used with touchless dispenser)

4 x 5 L canister

12 x 250mL Airless Spray

Also available in non-sterile packs:

12 x 0.5 L trigger spray

12 x 1 L capped bottle with a separate trigger

12 x 1 L bottle with flip cap

2 x 5 L canister