



ClearKlens

Tego® 2000 SS

VH25S

## Sterile ready to use amphoteric disinfectant with WFI

### Description

ClearKlens Tego® 2000 SS is a ready to use surface active amphoteric disinfectant for the pharmaceutical industry, presented in a 900 ml trigger spray.

### Key properties

- ClearKlens Tego® 2000 SS has been developed for the disinfection of walls, floors, work surfaces and process equipment within pharmaceutical cleanrooms.
- ClearKlens Tego® 2000 SS is a 900 ml ready to use solution presented in a bag in bottle delivery system to eliminate product contamination during use.
- ClearKlens Tego® 2000 SS is filtered to 0.2 µm and aseptically filled.
- ClearKlens Tego® 2000 SS is double bagged for transfer into the cleanroom.
- Tego® 2000 is the registered trademark of Evonik, Germany.

### Benefits

- Broad spectrum of activity
- Non tainting
- Odourless
- Low toxicity
- 0.2 µm filtered
- Double bagged for easy transfer into the cleanroom
- Sterility tested – Certificate provided

### Use instructions

ClearKlens Tego® 2000 SS is a ready to use 1% v/v aqueous solution with water for injection (WFI).

### Technical data

Appearance	Clear, colourless liquid
Relative density at 20°C	1
pH (neat solution at 20°C)	7.3
Chemical Oxygen Demand (COD)	4 gO <sub>2</sub> /kg
Biochemical Oxygen Demand (BOD)	0 gO <sub>2</sub> /kg
Nitrogen Content (N)	0 g/kg

*The above data is typical of normal production and should not be taken as a specification.*



Pharma™



ClearKlens

Tego<sup>®</sup> 2000 SS

VH25S

#### Safe handling and storage information

Store in original closed containers, away from sunlight and extremes of temperature. Full guidance on the handling and disposal of this product is provided in a separate Safety Data Sheet.

#### Product compatibility

ClearKlens Tego<sup>®</sup> 2000 SS when applied at the recommended concentration and temperature is suitable for use on the grades of stainless steel commonly found in the pharmaceutical industry. In the event of uncertainty it is advisable to evaluate individual materials before any prolonged use.

#### Test method

Validated titration test method available.

Each lot is provided with a Certificate of Analysis, available on [www.clearklens.com](http://www.clearklens.com)

#### Microbiological data

ClearKlens Tego<sup>®</sup> 2000 SS has been extensively tested using EN 1276, EN 1650 and EN 13697. It passes the EN14476 for MVA. These have shown ClearKlens<sup>®</sup> Tego 2000 SS to have a broad spectrum of activity against both Gram-positive and Gram-negative bacteria, yeasts and viruses.

#### Available packsize

6 x 900 mL trigger spray