

Product Information

Material	Polypropylene (PP) 
Construction	Nonwoven - meltblown
Saturant	70% IPA (USP Grade) with 30% DI water
Sterility	Gamma irradiated to a validated SAL 10 ⁻⁶
Shelf life	Sterile: 2 years from manufacturing date
Environment	ISO 5-8 Grade A/B for sterile, C/D for nonsterile

Efficacy Information

Test	Description	Log Reduction	Time	Test	Description	Log Reduction	Time
EN16615	<i>E. hirae</i>	>7.00	1 min	EN16615	<i>P. aeruginosa</i>	>7.19	1 min
EN16615	<i>S. aureus</i>	>6.62	1 min	EN16615	<i>C. albicans</i>	>6.36	3 min

Instructions for use

Contec PROSAT Meltblown Polypropylene Wipes with 70% IPA are ready-to-use.

When transferring the pouches to the point of use, remove each packaging layer as the environment becomes more critical.

Ensure the surface is uniformly covered with the solvent then wipe to dry with a Contec sterile cleanroom wipe.

Ordering Information

Description	Part No.	Size	Packaging
PROSAT Meltblown Polypropylene Wipes Presaturated with 70% IPA and deionised water	PS-850-BPR	200 x 200mm	50 wipes per pouch; 50 pouches per case
	PS-911-BPR	230 x 280mm	30 wipes per pouch; 50 pouches per case
	PSPP0043-BPR	230 x 280mm	50 wipes per pouch; 40 pouches per case
	PS-LPP-703-BPR	400 x 430mm (quarterfold)	50 wipes per pouch; 15 pouches per case
PROSAT Sterile Meltblown Polypropylene Wipes Presaturated with 70% IPA and deionised water	PS-840IR-BPR	200 x 200mm	40 wipes per pouch; 24 pouches per case
	PS-911EB-BPR	230 x 280mm	30 wipes per pouch; 48 pouches per case
	PSPP0039-BPR	230 x 280mm	50 wipes per pouch; 24 pouches per case
	PS-LPP-7030IR-BPR	400 x 430mm (quarterfold)	50 wipes per pouch; 12 pouches per case
PROSAT Sterile Meltblown Polypropylene LE Wipes Presaturated with 70% IPA and deionised water Less than 20EU per wipe	PS-911LE-BPR	230 x 280mm	30 wipes per pouch; 36 pouches per case

Use biocides safely. Always read the label and product information before use.

IAB Reinraum-Produkte GmbH
Erzberg 5
38126 Braunschweig
Deutschland

Telefon: +49 531 28484-0
Telefax: +49 531 28484-40
E-Mail: vertrieb@iab-bs.de